


High-Performance Rubber Seals Under Realistic Application Conditions with the Help of DMA

Xu Liang and Sascha Sebastian Riegler, Applications Laboratory Shanghai and Selb

Introduction

Acrylonitrile butadiene rubber (NBR, structural formula in figure 1) is a copolymer produced by the polymerization of acrylonitrile and butadiene monomers. The main process used to manufacture this caoutchouc is low-temperature emulsion polymerization [1]. The acrylonitrile content of the copolymers is typically between 18 and 50 mol.-% [1]. NBRs generally show good resistance to nonpolar solvents, high abrasion resistance, gas impermeability and good temperature resistance. As a result, they are widely used in the manufacture of various oil-resistant rubber products, such as bellows, gaskets and other seals, rubber gloves, oil-resistant soles, printing blankets, etc., and have become an indispensable elastic material in the automotive, aviation, petroleum, packaging, food, printing and other industries [2].

Some NBR products are subjected to constant strain and elevated temperatures during service. Therefore, knowledge of the relaxation and deformation set – either tension or compression set – is important to the customer during product design. When a material is used under constant strain, the material response may become irreversible under longer time scales and/or



1 Structural formula of acrylonitrile butadiene rubber [3].

higher temperatures. This can result in a non-zero, permanent deformation of the material after the strain is removed. This non-reversible part is an important factor in determining the applicability of certain rubber materials. There are several international standards and Chinese standards for testing the relevant relaxation and deformation set properties of elastomers, such as ASTM D395, GB/T 7759.1, GB/T 7759.2, and GB/T 1683.

However, information on the material performance for these properties can also be obtained with the NETZSCH DMA 303 *Eplexor*® by simulating the material behavior under application-relevant conditions.

APPLICATIONNOTE High-Performance Rubber Seals Under Realistic Application Conditions with the Help of DMA

Relaxation and Compression Set Measurements on As-Received and Post-Vulcanized NBR Experimental

Two different NBR samples were measured in compression mode with the DMA 303 *Eplexor*[®] using the appropriate compression steel sample holder and pushrod, as shown in figure 2. One is an as-received NBR sample that underwent a primary vulcanization process at 170°C in static air, and the other is a post-vulcanized NBR sample that was further heat-treated at 170°C for 2 h in an oven in static air. The diameter of the samples was 5.18 mm and 5.22 mm for the as-received and post-vulcanized NBR samples, respectively. The height of the sample was determined by the automatic length detection feature of the DMA 303 *Eplexor*[®].

The experiment was performed using the following six-segment procedure:

1. A static force of 0.05 N was applied to ensure contact with the sample during isothermal stabilization at 25°C for 5 min. At the end of the segment, the initial thickness, L_0 , was measured.
2. The temperature was then raised to 100°C at a heating rate of 10 K min⁻¹.
3. To stabilize temperature and allow the entire sample to equilibrate at 100°C, the temperature was held for 5 min before the next step.
4. A target static strain of 25% based on the length measured at the end of the previous segment was applied. The strain was held constant at this temperature for 60 min and the decay of the force and the relaxation modulus were observed as a function of time throughout the segment.
5. The applied force was reduced to the previous 0.05 N and then cooled back to 25°C at 10 Kmin⁻¹.
6. The temperature was held constant at 25°C for 20 min to stabilize the temperature and allow the sample to fully equilibrate at the given temperature. At the end of the segment, the length of the sample, L_1 , was measured again and the residual, non-reversible strain, $\epsilon = (L_1 - L_0)/L_0$, was determined.

The sample length measured at the end of the first segment is $L_0 = 7.722$ mm. After applying a static strain of -25 % at the beginning of the isothermal segment at 100°C, the static force decreases from its maximum value of 24.97 N to 20.41 N after one hour. Correspondingly, the relaxation modulus decreases from 4.77 MPa to 3.87 MPa. At the end of the measurement, the specimen has a length of $L_1 = 7.464$ mm. This corresponds to a residual strain of $\epsilon = -3.34$ % after one hour.

For the post-vulcanized NBR sample, a length of $L_0 = 7.638$ mm was measured prior to the start of the heating segment. The static strain of -25% requires an initial force of 21.41 N, which decreases to 17.10 N after 1 h at 100°C. The relaxation modulus decreases from an initial value of 4.06 MPa to 3.19 MPa during the isothermal segment. At the end of the experiment, a specimen length of $L_1 = 7.509$ mm was measured. Hence, the calculated residual strain in this case was $\epsilon = -1.69$ %.



2 Position of the sample between the sample holder and pushrod for compression-mode measurements


APPLICATIONNOTE High-Performance Rubber Seals Under Realistic Application Conditions with the Help of DMA

Measurement Results


While the as-received NBR still shows a residual strain of -3.34%, the post-vulcanized NBR sample only shows a value of -1.69%. This shows a drastic influence of the post-vulcanization process treatment on NBR highlighted by the reduction of the residual strain by about 50.6% compared to the as-received state.

From a microstructural point of view, the difference in residual strain can be explained by the higher degree of intermolecular chemical cross-linking of the polymer

chains for the post-vulcanized NBR sample. As a result, their mobility and ability to undergo configurational changes at elevated temperatures and/or under longer time scales is drastically reduced. Since irreversible, viscous flow requires movement of the main polymer chains into new metastable configurations, the increased degree of chemical cross-linking reduces the possibilities for configurational changes during deformation of the sample. The irreversible microstructural changes are reflected on the macroscopic scale by the reduction in force during the isothermal relaxation segment as shown in figures 3 and 4.



3 The diagrams show the experimental results of the NBR sample in the as-received state. They contain the values for the detected static force, F_{stat} , the average length, L_m , the relaxation modulus, E , and temperature as a function of time.



4 The diagrams show the experimental results for the NBR sample in the post-vulcanized state. They contain the values for the detected static force, F_{stat} , the average length, L_m , the relaxation modulus, E , and temperature as a function of time.

APPLICATIONNOTE High-Performance Rubber Seals Under Realistic Application Conditions with the Help of DMA

For product designers, the advantage of post-vulcanization of elastomers is that they can expect fewer physical and chemical changes within their product during service, such as the residual strain shown here. This allows them to tailor their final product more closely to the application of the material.

Conclusion

In addition, compared to relaxation and compression set experiments carried out in accordance with several international standards, dynamic mechanical analysis also allows in-situ observation of the reduction in force during constant strain. This can provide a product designer with additional information about the behavior of their material in service.

Literature

- [1] „Acrylnitril-Butadien-Kautschuk“. [Online]. Available under: <https://www.rado.de/spezialitaeten/nbr/>
- [2] „NBR: Acrylonitrile-butadiene rubber“. 17. Juli 2024. [Online]. Available under: <https://polymers.netzsch.com/Materials/Details/57>
- [3] „Acrylnitril-Butadien-Kautschuk“. 17. Juli 2024. [Online]. Available under: https://en.wikipedia.org/wiki/Nitrile_rubber