

# APPLICATION SHEET

Organics · Food  
STA 409 - SKIMMER




## Tobacco

### Introduction

Tobacco leaves are often smoked in form of a cigarette or cigar, or in a smoking or water pipe. Tobacco is also chewed and sniffed. Along with the pleasures of tobacco consumption come real risks of serious diseases such as lung cancer, respiratory and heart diseases. Major hazards of tobacco use involve carcinogenic compounds in tobacco

and tobacco smoke. One group of most probably dangerous substances are polycyclic aromatic hydrocarbons such as naphthalin, fluorine, anthracene, benzopyrene and benzoapyrene. For many people, it is also difficult to quit smoking. Many jurisdictions have thus enacted smoking bans in an effort to minimize possible damage to public health caused by tobacco smoking.



### Test Conditions

Temperature range: RT... 800°C  
Heating rate: 10 K/min  
Atmosphere: Argon at 75 ml/min  
Sample mass: 578 mg  
Crucible: Al<sub>2</sub>O<sub>3</sub> beaker  
Sensor: TGA type S

### Test Results

Two tobacco samples were characterized by thermogravimetry (full lines) and mass spectroscopy (dashed lines). For clarity, only particular mass number 252 is shown which is most probably due to benzoapyrene. It is, however, important to emphasize that the mass spectrum of tobacco contains a large amount of mass numbers which are due to tar, nicotine derivatives and other compounds mentioned above. The measurements shown in the diagram revealed clear differences between the two tobacco samples: sample "Control" exhibited a larger mass loss at low temperatures as well as a stronger evolution of the polycyclic aromatic hydrocarbon benzoapyrene.