

APPLICATION SHEET


Inorganics · Chemical Industry
TG 209 F1 Libra®

Aluminum Hydroxide

Introduction

Aluminum hydroxides are the main compounds of bauxite. The principal aluminum hydroxide minerals found in varying proportions with bauxites are gibbsite ($\text{Al}(\text{OH})_3$) and the polymorphs boehmite ($\alpha\text{-AlOOH}$) and diaspore ($\beta\text{-AlOOH}$). The bulk of world bauxite production (approximately 85%) is processed into aluminum oxide (Al_2O_3)

via a wet chemical, caustic leach method (Bayer process). The resulting Al_2O_3 is then reduced to aluminum metal (Al) using an electrolytic process, the Hall-Heroult processor used in high-temperature ceramic applications. Purified aluminum hydroxides and oxy-hydroxides were used as flame-retardent fillers in polymers, employed as paper coating pigments in paints and varnish.



Test Conditions

Temperature range: RT ... 1000°C
Heating rate: 10 K/min
Atmosphere: Air at 40 ml/min
Sample mass: 36 mg
Crucible: Alumina
Sensor: Platinel

Test Results

The two $\text{Al}(\text{OH})_3$ samples are nearly pure indicated by the mass loss of 34.7% and 34.4%. The theoretical mass loss from $\text{Al}(\text{OH})_3$ to Al_2O_3 is 34.6%. This means there might be surface water at sample 1 and some boehmite amounts with the second one. Sample 1 additionally shows a three-step degradation in comparison with sample 2.