

APPLICATION SHEET


Inorganics · Building Materials
ST 449 **F1 Jupiter**[®]

Cement / Gypsum

Introduction

Portland cement is the most common type of cement. It consists of a mixture of oxides of calcium, silicon and aluminium. The main phases are tricalcium silicate (C_3S), dicalcium silicate (C_2S), tricalcium aluminate (C_3A) and a ferrite

phase C_4AF . Portland cement is produced by heating limestone with clay and/or sand up to about 1480°C . To the resulting clinker, 4-5% gypsum is mixed and then grinded and milled to a fine grain size. The additives gypsum, anhydrite etc. influence the setting time of the cement.



Test Conditions

Temperature range: RT ... 550°C
Heating rate: 10 K/min
Atmosphere: Air (50 ml/min)
Sample mass: 39.5 mg
Crucible: Al with 50 μm hole
Sensor: TGA-DSC type S

Test Results

The separation and quantification of calcium sulfate dihydrate (DH) and hemihydrate (HH) in cement samples is only possible when a certain water vapor pressure over the sample can be created. A possibility is to use closed aluminum crucibles with a small pinhole of about 40 to 50 μm . From the TGA curve, the DH amount can be calculated to 2.5% and the HH amount to 1.8%. The TGA step (0.5%) starting at about 420°C is due to the decomposition of $\text{Ca}(\text{OH})_2$ (content approx. 2.2%).