

# APPLICATION SHEET


Polymers · Polymer Manufacturing  
TG 209 **F1 Iris**® – FT-IR

## Polyethylene

### Introduction

Polyethylene (PE) is a thermoplastic material which is heavily used for consumer products as foils (wrapping,

packaging), containers (bottles, tanks), pipes, tubes or other engineered products. PE is odorless, flavorless, physiological indifferent and can therefore be used in the food industry.




### Test Conditions:

Temperature range: RT ... -165°C  
Heating/cooling rates: 20 K/min  
Atmosphere: Nitrogen at 40 ml/min

Sample mass: 7.7 mg  
Crucible: Alumina  
Sensor: Platinel

# APPLICATION SHEET

Polymers · Polymer Manufacturing  
TG 209 **F1 Iris**® – FT-IR



## Test Results

Figure 1 depicts the 3-dimensional plot of the FT-IR spectra and the TGA curve. PE evolves completely, starting at

about 420°C. The pyrolysis products are hydrocarbon fragments showing C-H and C-C bondings in the FT-IR (figure 2)