

## Binder Burnout of Stabilized Zirconia

Dr. Michael Schöneich and Senol Gezgin

## Introduction


Zirconium dioxide (ZrO<sub>2</sub>) is one of the most often studied ceramic materials. Upon heating, zirconia undergoes disruptive phase changes. By adding small percentages of yttria, these phase changes are eliminated, and the resulting material has superior thermal, mechanical and electrical properties.

## **Measurement Conditions and Results**

This STA measurement between room temperature and 1200°C exhibits two small losses up to 450°C (3.4% in total; green TGA curve) which correspond very well with the two exothermic peaks at 197°C and 399°C in the DSC curve (blue). These effects (mass loss up to 500°C, exothermic peaks with high enthalpies) can be typically observed during the binder burnout of ceramic materials. The small endothermic DSC peak at around 67°C is caused by the melting of the binder.



1 STA 449 *F5* Jupiter®



2 STA measurement of zirconium dioxide (26.2 mg) in platinum crucibles at a heating rate of 10 K/min

