


## Polyetherimide

## Introduction

Polyetherimide is an aromatic, amorphous thermoplastic containing both ether links and imide groups. It has a good fire resistance and a good thermal stability. Furthermore, it

is chemical resistant and features high dielectric strength and extremely low smoke generation. Is is often used for heat-resistant products such as in microwave ovens and circuit boards.



## **Test Conditions**

Temperature range: -25 ... 400°C Heating rate: 10 K/min

Atmosphere: Nitrogen (20 ml/min)

Sample mass: 9.37 mg
Crucible: Al, pierced lid

## **Test Results**

The endothermic step at 225.2°C (midpoint) with a change in specific heat of 0.21 J/(g·K) is due to the glass transition of polyetherimide. The increase in the specific heat flow rate versus temperature is due to the increase in specific heat.

