APPLICATION NOTE

TGA-FT-IR Coupling – Improved Sensitivity by Means of Fast Heating Rates

Dr. Ekkehard Füglein and Sonja Eichholz

Introduction

The field of thermal analysis comprises methods for the characterization of physical and chemical properties or property changes as a function of temperature. Thermogravimetry allows for the quantification of mass changes, for example, the release of reaction and decomposition gases. When these gases are transferred into a gas measuring cell, identification of the released gases becomes also possible. The so-called TGA-FT-IR coupling is a tried-and-tested combination of an analytical and a spectroscopic analysis method.

A quite new addition to the established STA 449 **F1** *Jupiter*[®] (figure 1) is the high-speed furnace (cross-section in figure 2) which can be operated at heating rates of up to 1000 K/min (furnace systems for the most varied of applications covering a temperature range from -150°C to 2400°C are now available).



2 Cross-section of the high-speed furnace

The influence of the heating rate and the related release rate on the thermogravimetric and spectroscopic measurement results will be discussed in this application note.



1 STA 449 F1 Jupiter[®] with Tensor 27[™]



APPLICATIONNOTE TGA-FT-IR Coupling – Improved Sensitivity by Means of Fast Heating Rates

Results

a) Polypropylene PP

When the heating rate during thermoanalytical experiments is varied, the detected effects are shifted to higher temperatures with increasing heating rates (figure 3). This is well-known and can be used for kinetic evaluations. With an increasing release temperature, the release rate

also increases significantly (figure 4). Therefore, the concentration of the sample gases to be analyzed in the constant carrier gas flow also rises and the sample gases can easily be detected and identified. The mass-loss steps, however, are not dependent on the heating rate.



3 TGA results for polypropylene (PP)



4 DTG results for polypropylene (PP)



APPLICATIONNOTE TGA-FT-IR Coupling – Improved Sensitivity by Means of Fast Heating Rates

b) CaCO₃

The relation between the heating rate and decomposition temperature at a constant step height, discussed for the pyrolysis of propylene, can also be seen during the thermal decomposition of calcium carbonate into calcium oxide and carbon dioxide (figures 5 and 6).



5 TGA results for calcium carbonate (CaCO₃)



6 DTG results for calcium carbonate (CaCO₃)



APPLICATIONNOTE TGA-FT-IR Coupling – Improved Sensitivity by Means of Fast Heating

Rates

Figure 7 shows the absorption intensity of the corresponding Gram-Schmidt traces which is expected to increase with increasing heating rates. It should be noted here that the transport of the released sample gases into the IR gas measuring cell is hardly delayed due to the fast heating rates. This can be seen by the comparison of the maximum release rate (DTA) with the maximum IR intensity (GS) in figure 8.



7 Gram-Schmidt results, calcium carbonate (CaCO₃)



8 Comparison of the DTG temperature with the IR intensity (GS)



APPLICATIONNOTE TGA-FT-IR Coupling – Improved Sensitivity by Means of Fast Heating Rates

c) $CaC_2O_4 \times H_2O$ mixed with SiO₂

For investigation of the detection limit, a mixture of calcium oxalate monohydrate (CaC₂O₄ x H₂O) and quartz sand (SiO₂) was prepared. The selected mixing ratio was 1:10 so that the expected release of water corresponded to approximately 1% of the sample mass.

The thermal release of approximately 1% of water from this mixture could not be detected at a heating rate of 20 K/min; employing a heating rate of 200 K/min, however, it could clearly be detected (figure 9 to 11).



TGA results and Gram-Schmidt trace (dasehd line) of a sample of calcium oxalate monohydrate and SiO, as 9 inert matrix



10 Presentation of all IR spectra at a heating rate of 200 K/min



APPLICATIONNOTE TGA-FT-IR Coupling – Improved Sensitivity by Means of Fast Heating Rates



11 Enlarged scaling of figure 11 with the absorption bands for water

Summary

By means of fast heating rates of up to 500 K/min, it is possible to considerably increase the release rates of gasesous products from a sample. Therefore, their concentration also increases in comparison with the carrier gas, resulting in a significant improvement of the detection limit of a TGA-FT-IR coupling.

