

About the Thermal Behavior of Acetylsalicylic Acid and Aspirin®

Dr. Ekkehard Füglein¹, Dr. Carolin Fischer¹ and Dr. Albrecht Rager²

¹NETZSCH-Gerätebau GmbH, Wittelsbacherstr. 42, 95100 Selb

²Bruker Optik GmbH, Rudolf-Plank-Straße 27, 76275 Ettlingen

Introduction

Early high cultures already used willow bark as a remedy for fever and pain [1]; the Roman scholar, Plinius the Elder, regarded willow bark as a medicine and the Teutons and Celts produced extracts by cooking willow bark, the ingredients of which were chemically related to synthetic acetylsalicylic acid [2]. Although various chemists were able to produce salicin and salicylic acid in the 19th century, it was not until 1897 that Felix Hoffmann succeeded in synthesizing acetylsalicylic acid without impurities at BAYER's headquarters in Wuppertal-Elberfeld, Germany. Kurt Wittauer (figure 2) tested this drug on patients in the following years until BAYER (figure 1) finally filed for the corresponding patent in 1921. The painkiller began its triumphant success around the world and today, BAYER produces more than 50,000 tons of acetylsalicylic acid per year [4].



1 BAYER Aspirin®


Aspirin, a New Salicyl Preparation

By Dr. Kurt Wittauer, Senior Physician at the Deaconesses' House Halle a. S.

Nowadays it takes a certain amount of courage to recommend a new remedy. They are thrown onto the market almost every day, and you would have to have a great memory if you wanted to keep all the new names in mind. Many appear, are praised and recommended by individual authors and especially by the factories, and after a short time, you don't hear of them anymore.

When, almost a year ago, the Farbenfabriken vorm. Friedr. Bayer & Co. in Elberfeld, Germany, sent me a new salicyl drug for testing, I approached its application with a certain degree of mistrust.

The new drug, called Aspirin, a compound of acetic and salicylic acid, has the chemical formula



and forms white crystal needles with a melting point of 135°C, which dissolve in water at 37°C to 1%. It is easily soluble in typical organic solvents; with iron chloride, the solution no longer appears blue in color.

It dissolves easily in diluted alkalis and breaks down there into its two components after a few minutes, so the conclusion that the substance is only split into its constituents when in contact with alkaline liquids – i.e., in the intestine – seems justified. Of course, where there is a lack of hydrochloric acid in the stomach, as would be the case with carcinoma present, the decomposition will take place already in the gastric juice.

The most important difference between aspirin and ordinary salicylic acid would be that it does not attack the stomach, passes through it unchanged and only cleaves in the alkaline intestinal juice.

Prof. Dreser will be giving pharmacological information about the drug elsewhere.

2 Original document by Kurt Wittauer [3] in German (1899), translated into English by D. Steidl/N. Huss for this Application Note

Drugs containing the active ingredient acetylsalicylic acid are available in various pharmaceutical forms and are employed not only because of their analgesic effect but also due to their anti-inflammatory, antipyretic and antiplatelet properties.

Pure acetylsalicylic acid is a pure white powder that is poorly soluble in water, has a melting point of 136°C and decomposes at higher temperatures. Various methods of thermal analysis, infrared spectroscopy, and combinations of the two were employed in this work to investigate gaseous decomposition products.

APPLICATIONNOTE About the Thermal Behavior of Acetylsalicylic Acid and Aspirin®

Methods and Preparation

Acetylsalicylic acid (CAS: 50-78-2) was acquired from Sigma Aldrich with a purity of > 99%. For the investigation of the original substance, the BRUKER TENSOR II was used to measure the samples with attenuated total reflection (ATR). For determination of the melting

behavior, the NETZSCH DSC 214 *Polyma* was used. For thermal characterization of the gases released, a thermobalance was coupled to an infrared spectrometer – the NETZSCH TG 209 *F1 Libra*® to the Bruker Equinox 55/S. The measurement conditions for the thermo-analytical and spectroscopic investigations are summarized in tables 1 to 3.

Tab 1. Measurement conditions for the DSC investigation of acetylsalicylic acid

Acetylsalicylic Acid	
Sample mass	2.08 mg
Crucible material	Aluminum, pierced
Crucible mass	52.75 mg
Temperature range	25 ... 160°C
Heating rate	7 K/min
Atmosphere	Nitrogen (50 ml)

Tab 2. Measurement conditions for the thermogravimetric investigation of an Aspirin® tablet by means of TGA-FT-IR

Aspirin®	
Sample mass	9.141 mg
Crucible material	Alumina, open
Crucible mass	162.75 mg
Temperature range	25 ... 600°C
Heating rate	10 K/min
Atmosphere	Nitrogen (40 ml)
Scans	32
Resolution	4 cm ⁻¹
Spectral range	650 - 4500 cm ⁻¹


Tab 3. Measurement conditions for the spectroscopic investigation of (ATR) of acetylsalicylic acid

Acetylsalicylic Acid	
Detector	DTGS
Scans	32
Resolution	4 cm ⁻¹
Spectral range	650 - 4500 cm ⁻¹


Results and Discussion

Investigation of the active ingredient acetylsalicylic acid with the help of FT-IR spectroscopy yields an infra-red spectrum at room temperature which is in good agreement with the library spectrum (Bruker ATR-LIB-Polymers-1-472-2) (figure 3). The melting range of acetylsalicylic acid is indicated by the manufacturer to be at 134°C to 136°C. The investigation by means of Differential Scanning Calorimetry (DSC) provides a melting

enthalpy of 178 J/g and a temperature for the extrapolated onset of 138.5°C. As can also be clearly seen from figure 4, the heat-flow signal indicates the beginning of the melting process of the sample already at significantly lower temperatures than determined by the standard-compliant evaluation for the extrapolated onset. In literature, two polymorphic forms of acetylsalicylic acid are described: Form I with a melting temperature of 144.9°C and Form II with a melting temperature of 135.5°C [5, 6].



3 FT-IR spectrum of acetylsalicylic acid (red) by means of attenuated total reflectance (ATR) and comparison with a library spectrum (Bruker ATR-LIB-Polymers-1-472-2) of o-acetoxybenzoic acid (acetylsalicylic acid, blue)




4 DSC results for acetylsalicylic acid

APPLICATIONNOTE About the Thermal Behavior of Acetylsalicylic Acid and Aspirin®


Above about 150°C, the thermal degradation of acetylsalicylic acid begins. Therefore, thermogravimetry (TGA) is better suited for further characterization above the melting point (figure 5).

For characterization of the thermal degradation, a piece of an aspirin tablet was investigated with the help of the TGA-FT-IR coupling. Although the thermogravimetric results between 150°C and 450°C show a two-step thermal degradation reaction and the amounts of gases

released can be quantified, it is not possible to determine which gases are responsible for the detected mass loss without spectroscopic analysis. If one carries out a measurement where the thermobalance is coupled to an infrared spectrometer, the gas phase can be investigated continuously during the entire measurement. All infrared spectra are presented in a three-dimensional arrangement, temperature-scaled, in figure 6. The results of the thermogravimetric measurement can also be seen in the left rear area.



5 Thermogravimetric results of an Aspirin® tablet




6 Three-dreidimensional temperature-dependent presentation of all spectra of the measurement on the Aspirin® tablet with the TGA curve (red) in the rear area of the cube

APPLICATIONNOTE About the Thermal Behavior of Acetylsalicylic Acid and Aspirin®


If individual spectra are extracted from this presentation at the temperatures with the highest absorption intensities, the gases released can be identified with the help of comparison spectra from gas phase libraries. The individual spectrum for the first mass-loss step at 180°C, which is characteristic, is in very good agreement with the spectrum for acetic acid from the EPA-NIST gas phase library (figure 7). The red arrows indicate absorption bands that do not match acetic acid, but correspond very well with the absorption bands for salicylic acid (EPA-NIST). This leads to the assumption that the acetylsalicylic acid, as in reaction equation 1, thermally degrades into salicylic acid and acetic acid (Equation 1). At 180°C, the acetic acid formed is already gaseous while the salicylic acid, with a melting point of 159°C, begins to evaporate. This is certainly also the reason why the first mass-loss step passes directly into the following step. The combination of decomposition and evaporation confirms the degradation mechanism proposed by Rebeiro et al. [7].

In conjunction with the tablet form of the active ingredient acetylsalicylic acid, the influence of humidity on the reaction products of the thermal degradation is emphasized along with additives such as starch and magnesium stearate monohydrate. Gupchup et al., however, point out that the dry active ingredient acetylsalicylic acid can itself ensure the presence of water through dimerization in the sense of condensation [8].

When comparing the two spectra for acetic acid and salicylic acid, it is noticeable that the absorption bands in the range between 1760 cm⁻¹ and 1820 cm⁻¹ can only be attributed to acetic acid, while the absorption bands between 1460 cm⁻¹ and 1500 cm⁻¹ represent salicylic acid. If the intensity course of the absorption ranges is calculated as a function of temperature, "traces" are obtained for each substance; these are proportional to the corresponding amounts released as a function of temperature.



7 Extracted individual spectrum of an Aspirin® tablet at 180°C (red) and comparison with the spectra from the gas phase library (EPA-NIST) for acetic acid (black) and salicylic acid (blue)




APPLICATIONNOTE About the Thermal Behavior of Acetylsalicylic Acid and Aspirin®


A comparison of these temperature-dependent traces for acetic acid and salicylic acid is shown in figure 8 with the Gram-Schmidt trace (sum of the intensities that are not dependent upon wavelength) and the TGA signal. As with the TGA signal, the Gram-Schmidt trace reveals that the first mass-loss step passes directly and without plateau into the second mass-loss step. The reason for this can be found in the traces of the two products, which show that the release of acetic acid can be detected up

to about 300°C and, in addition, the evaporation of salicylic acid starts already at lower temperatures.


Along with salicylic acid, also the formation of carbon dioxide can be detected by means of the temperature-dependent course of the absorption intensities. This is confirmed by the extracted individual spectrum at 360°C (figure 9).



- 8 TGA results for an Aspirin® tablet (black) with the Gram-Schmidt trace (blue) as well as the trace for salicylic acid (green dashed), acetic acid (red dashed) and carbon dioxide (black dashed)



- 9 Extracted individual spectrum of an Aspirin® tablet at 360°C (red) and comparison with the library spectrum of salicylic acid (green)



In the range between wave numbers 2424 and 2224, the absorption bands of CO₂ are clearly visible. In addition, there are indications that phenol has formed. The positions of the most intense absorption bands for phenol are marked with red arrows. It can therefore be assumed that – along with the evaporation of salicylic acid – a decomposition process also takes place; this suggests the formation of phenol and CO₂ as shown in equation 2.

Summary

Acetylsalicylic acid was investigated using FT-IR spectroscopy at room temperature (ATR), and the FT-IR spectra obtained were used for identification by means of comparison with a spectra library. The DSC was used for investigation of the melting behavior. Additionally, the thermal behavior of Aspirin® was characterized by means of TGA-FT-IR. The spectra for the gases released during thermal treatment were compared with a gas phase library for identification of the products. It was thus possible to confirm degradation mechanisms known from literature and it was additionally shown that the common additives used in the tableting of Aspirin® appear to have no detectable influence on the formation of gaseous decomposition products.

Literature

- [1] K.-C. Bergmann, J. Ring (Hrsg.), *History of Allergy*. Verlag Karger Medical and Scientific Publishers, 2014, S. 132
- [2] O. Lafont, "From the willow to aspirin", *Rev Hist Pharm. (Paris)*. **55** (354), Jul 2007, S. 209–216.
- [3] K. Wittauer, „Aspirin, ein neues Salicylpräparat“, *Therapeutische Monatshefte* **13**, S. 330 (1899)
- [4] BAYER, Chemie der Medikamente, Unterrichtsmaterialien Galenik, 4a
- [5] P. Peksa, J. Trzmiel, M. Ptak, M. Kostrzewa, R. Szatanik, A. Barascu, D. Enke, A. Sieradzki, „Confinement-induced polymorphism in acetylsalicylic acid–nanoporous glass composites“, *J. Mater Sci*, **54** (2019) 404 doi.org/10.1007/s10853-018-2853-8
- [6] A.G. Shtukenberg, C.T. Hu, Q. Zhu, M.U. Schmidt, W. Xu, M. Tan, B. Kahr, „The Third Ambient Aspirin Polymorph“, *Cryst. Growth Des.*, **17** (2017) 3562–3566, DOI: 10.1021/acs.cgd.7b00673
- [7] Y.A. Ribeiro, A.C.F. Caires, N. Boralle, M. Ionashiro, „Thermal decomposition of acetylsalicylic acid (aspirin)“ *Thermochimica Acta* **279** (1996) 177-181.
- [8] G. Gupchup, K. Alexander, D. Dollimore, „The use of thermal analysis and mass spectrometry to study the solid state behavior in pharmaceutical tablet mixtures“ *Thermochimica Acta* **196** (1992) 267-278