# APPLICATIONNOTE

# Challenging Samples: A Solution to Monitor Curing Use of the Kinexus Disposable Plate Systems

Dr. Adrian Hill

# Introduction

The Kinexus rotational rheometer system by NETZSCH has a range of disposable plate systems (Table 1) all based on the same common design (see Figure 1), for the various different temperature control systems available. These are perfect for when the sample is expected to cure/adhere strongly to the measuring system, making standard cleaning all but impossible. For these challenging samples, the unique Kinexus design enables routine measurements with minimal material waste/ environmental impact and low cost with a rapid test turnaround time.

This system comprises of a new quick connect (non-disposable) upper geometry shaft (Figure 1, (1)) with autorecognition/configuration and a lower plate (2), which

are both used to securely hold the low-priced and disposable sample facing plates (4). The upper plates (3) are available in a variety of sizes (typically from 10 mm to 40 mm, see Table 1) in various quantities depending on testing needs.

# **Kinexus Active Hood**

The Kinexus active hood contains unique technology to **prevent** the formation of thermal gradients across the sample. For measurements where the **absolute** temperature is key, the Kinexus active hood system should be considered.



1 Kinexus disposable plate system (KNX2155)



Use of the Kinexus Disposable Plate Systems

Table 1 and Table 2 show the different cartridge systems and disposable plates available for the Kinexus rotational rheometers.

# rSpace for Kinexus: Rheology Software to Guide through the Measurements

rSpace takes a unique guided approach to rheology measurements with test methods are called sequences.

#### Table 1. Disposable plate system compatibility

Kinexus Cartridge System	Temperature Range
Plate (Peltier)	-40 (-5)* to 200°C
Plate (Peltier), -E variant	-40 (-5)* to 150°C
Active hood (Peltier)	-40 (-5)* to 200°C
Active hood (Peltier), -E variant	-40 (-5)* to 150°C
HTC plate	0 (ambient)* to 350°C
Active hood (Peltier) Active hood (Peltier), -E variant HTC plate	-40 (-5)* to 200°C -40 (-5)* to 150°C 0 (ambient)* to 350°C

\* depending on the otional cooling system

 Table 2.
 Common disposable plate consumables

Disposable Plate Consumables	Material
Disposable plate sets (upper and lower)	
40-mm set	Aluminum upper, stainless steel lower
40-mm set	Stainless steel upper and lower
25-mm set	Aluminum upper, stainless steel lower
25-mm set	Stainless steel upper and lower
Disposable lower plate only	
Disposable lower plate	Stainless steel
Disposable upper plates only	
40 mm	Aluminum
25 mm	Aluminum
12 mm	Aluminum
10 mm	Aluminum

\* For details, please see the Kinexus accessory brochures

File Edit View Bun Sequence Jools Window r@ges Help       Favourites Favourites Results Favourites Favourites Favourites Results Favourites Favourites Favourites Favourites Favourites Results Favourites <	rSpace		
<ul> <li>• Inport_0003 Import Rosand shear and extensional data in</li> </ul>	<u>F</u> ile <u>Edit View R</u> un Sequence <u>T</u> ools <u>W</u> indow r	rPages Help	
Favourites       Import_0003 Single Frequency (Strain control)         Tookkt_0003 Single Frequency (Strain control)         Tookkt_0013 Single Frequency (Strain Control) <tr< td=""><td>🗋 🗣 📂 🔊 rFinder 🚽 + 📑 🗋 🖉 Chart 🔹 🛅 Table -</td><td>🝷 🗧 Report 👻 🧒 🗍 Zero gap 🗃 Load sample 🝷 🗃 Unload sample \mid 🔘 Start currer</td></tr<>	🗋 🗣 📂 🔊 rFinder 🚽 + 📑 🗋 🖉 Chart 🔹 🛅 Table -	🝷 🗧 Report 👻 🧒 🗍 Zero gap 🗃 Load sample 🝷 🗃 Unload sample \mid 🔘 Start currer	
Favourites       ×         Favourites       ×         Favourites       ×         * Tookit_V001 Shear Rate Table       *         * Tookit_V003 Yield Stress (Stress Ramp)       *         * Tookit_V004 Single Shear Rate       *         * Tookit_V005 Shear Rate Ramp       *         * Tookit_V005 Shear Rate Ramp       *         * Tookit_V005 Shear Rate Ramp       *         * Tookit_0001 Amplitude Table (Strain control)       *         * Tookit_0002 Frequency (Strain control)       *         * Tookit_0003 Single Frequency (Strain control)       *         * Tookit_0003 Single Frequency (Strain control)       *         * Tookit_0001 Creep and Recovery       *         * Tookit_R001 Stress Relaxation       *         * Tookit_R003 Stress Relaxation       *         * Tookit_R003 Stress Relaxation       *         * Tookit_R003 Stress Relaxa	🗧 🔶 Favourites 🕒 Results 📑 Go to top stop 📑 Sample de	details 🕂 Manual gap 👔 Set temperature 🔚 Viscometry 🔹 🚾 Oscillation 🕞 🖂 Creep &	
Favourites       Results            • Add to Favourites          • • • • • • • • • • • • • • •	Favourites X	a mater	
Add to Favourites X O A Add to Favourites X O A Add to Favourites X O Add to Favourite X O	Favourites Results	Za rrinder	
Tookkt_V001 Shear Rate Table         Tookkt_V002 Three Step Shear Rate         Tookkt_V003 Yield Stress (Stress Ramp)         Tookkt_V004 Single Shear Rate         Tookkt_V005 Shear Rate Ramp         Tookkt_0001 Amplitude Table (Strain control)         Tookkt_0003 Single Frequency (Strain control)         Tookkt_0004 Single Frequency (Strain control)         Tookkt_0001 Creep and Recovery         Tookkt_R001 Stress Relaxation         Timport_0001 Import data from Bohlin software         Import_0002 Import microrheology data in to rSpace         Import_0003 Import Rosand shear and extensional data in	🐤 Add to Favourites 🗙 😻 🔞 🕇 🦊		
Search for Crookit_V004 Single Shear Rate Sockit_V005 Shear Rate Ramp Tookit_V005 Shear Rate Ramp Tookit_0001 Amplitude Table (Strain control) Tookit_0002 Frequency (Table (Strain control) Tookit_0003 Single Frequency (Strain control) Tookit_0001 Creep and Recovery Tookit_R001 Stress Rehxation Tookit_R001 Stress Rehxation Tookit_0002 Import data from Bohlin software Timport_0003 Import Rosand shear and extensional data in	Toolkit_V001 Shear Rate Table	1 L 🕂 🚵 🥦 🖉 🖽 🔛 .	
Search for       disposable         Seample_0026 Disp	t Toolkit_V002 Three Step Shear Rate		
Incoded_v004 Single Single Single Kate         Tookit_v005 Shara Rate Ramp         Tookit_v005 Shara Rate Ramp         Tookit_0002 Amplitude Table (Strain control)         Tookit_0003 Single Frequency Table (Strain control)         Tookit_0004 Single Frequency (Strain control)         Tookit_0003 Single Frequency (Strain control)         Tookit_0004 Single Frequency (Strain control)         Tookit_0004 Single Frequency Temperature Ramp(Strain         Tookit_0001 Creep and Recovery         Tookit_0001 Import data from Bohlin software         Import_0002 Import microrheology data in to rSpace         Import_0003 Import Rosand shear and extensional data in	at Toolkit_V003 Yield Stress (Stress Ramp)	Search for Visposable	
Tookit_0001 Amplitude Table (Strain control)         Tookit_0001 Amplitude Table (Strain control)         Tookit_0002 Frequency Table (Strain control)         Tookit_0003 Single Frequency (Strain control)         Tookit_0004 Single Frequency (Strain control)         Tookit_0004 Single Frequency Temperature Ramp(Strain         Tookit_0011 Amplitude Table (Strain control)         Tookit_0004 Single Frequency (Strain control)         Tookit_0014 Single Frequency (Strain control)         Tookit_0014 Single Frequency Temperature Ramp(Strain         Tookit_0011 Stress Relaxation         Timport_0002 Import data from Bohlin software         Timport_0002 Import microrheology data in to rSpace         Tookit_0003 Import Rosand shear and extensional data in	Tookit_V004 Single Snear Rate	Location: To Date: To Date:	
Tookit_0002 Frequency Table (Strain control)         Tookit_0003 Single Frequency (Strain control)         Tookit_0004 Single Frequency Temperature Ramp(Strain         Tookit_001 Stress Relaxation         Tookit_0002 Import data from Bohlin software         Import_0002 Import microrheology data in to rSpace         Import_0003 Import Rosand shear and extensional data in	Toolkit 0001 Amplitude Table (Strain control)	rSpace v 01/01/2008 v 01/08/202	
**       Tookkt_O003 Single Frequency (Strain control)         **       Tookkt_O004 Single Frequency Temperature Ramp(Strain         **       Tookkt_O001 Creep and Recovery         **       Tookt, R001 Stress Relaxation         **       Tookt O002 Import data from Bohin software         **       Import_0002 Import microfheology data in to rSpace         **       Topolity Constant Stress Relaxation         **       Import_0003 Import Rosand shear and extensional data in	Toolkit_0002 Frequency Table (Strain control)		
Sample_0026 Disposable plate sample loading.rseq Sample_0026 Disposable plate sample loading.rseq Sample_0027 Disposable plate unload sample.rseq Tookt_C001 Import data from Bohlin software Tookt_0002 Import_0002 Import microfheology data in to rSpace Tookt_0003 Import_0003 Import Rosand shear and extensional data in	Toolkit_0003 Single Frequency (Strain control)	Name	
*** Tookkt_C001 Creep and Recovery         *** Tookkt_R001 Stress Relaxation         *** Tmoort_0001 Import data from Bohlm software         *** Import_0002 Import microrheology data in to rSpace         *** Import_0003 Import Rosand shear and extensional data in    Sample_0227 Disposable plate unload sample.rseq rSolution_0030 Disposable plate single frequency timed.rseq rSolution_0031 Disposable plate single frequency with temperature ramp.rseq rSolution_0031 Disposable plate single frequency blate	Sample_0026 Disposable plate sample loading.rseq           "st Toolkit_C001 Creep and Recovery           "st Toolkit_C001 Creep and Recovery           "st Toolkit_C001 Creep and Recovery		
rsolution_0030 Disposable plate single frequency timed.rseq rsolution_0031 Disposable plate single frequency with temperature ramp.rseq			
rSolution_0031 Disposable plate single frequency with temperature ramp.rseq	Import 0001 Import data from Boblin software	01 Stress Relaxation 01 Import data from Bohin software	
Import_0003 Import Rosand shear and extensional data in	Import_0002 Import microrheology data in to rSpace	rSolution 0031 Disposable plate single frequency with temperature ramp rseq	
	Import_0003 Import Rosand shear and extensional data in		

2 Included as standard, disposable plate sequences for guided use



Use of the Kinexus Disposable Plate Systems

Hundreds of pre-configured sequences are supplied with rSpace for users to use and edit for their own specific needs. This includes sequences for curing samples that are chemically\* initiated (e.g., two-part epxoy systems) or thermally\*\* initiated (e.g., thermosets). A custom "Event timer" is included so that samples that are initiated by an external means (for example mixing a two part-epoxy) are all timed (for accurate comparisons of different samples) to the critical time when the sample is mixed (see Figure 3). All disposable plate sequences guide the user step-bystep on how to assembly, mount, use and remove the disposable plate systems (Figure 4). Of course, this is all with the standard rSpace flexibility of being able to be customized for every specific testing requirement.

Design sequences, labelled "Sample\_" are also included to both separately load (Sample\_0026) and unload (Sample\_0027) the disposable plates and curing sample, or for use to be imported into your own custom measurement sequence.



# Usage Tip Depending on the particular sample properties, it is sometimes possible to remove (peel off) the flexible lower plate system along with the sample to enable

to reuse of the disposable upper plate.

3 Included control of the customer "Event time" to time when the reaction starts



4 Kinexus sequences guide through the use of the disposable plates

\* in Kinexus rSpace software rSolution\_0030

\*\* in Kinexus rSpace software rSolution\_0301





Use of the Kinexus Disposable Plate Systems

# Usage Tip

For samples that have a large tendency to shrink, it is advisable to place a drop of liquid (water for lower temperatures, oil for higher temperatures) between the lower plate adapter and the disposable lower plate. This further increases the rigidity of the system.

# Curing Profiles Measured by Rheometry

Typical curing profiles of two different samples measured with the disposable plate system are shown in Figure 5 using a single frequency oscillation. Typical properties are that the phase angle (green curve) starts high (i.e., viscous/liquid like), and then lowers (i.e., more elastic/ like) as the cure proceeds. Simultaneously, the materials (complex) modulus (as indicated by the combination of elastic and viscous modulus) increases as the sample gets "stiffer".

This can also be monitored with the elastic (G', red symbols) and viscous (G'', blue symbols) modulus where

## **Phase Angle**

The phase angle is a relative measure of the viscous and elastic properties of a material. It range from 0° for a fully elastic meaterial to 90° for a fully viscous material. It can be considered the "degree of fluidity".

they both tend to increase as the cure progresses, with initially the viscous modulus dominating (as more liquid like, phase angle >45°) then a "cross-over point" where the viscous modulus = elastic modulus (due to the phase angle = 45°) which can also be referred to as a "gel point"<sup>1</sup>, until finally the elastic modulus dominates over the viscous modulus. Depending on the material, a plateau in the curve of the phase angle and/or the moduli may indicate the final cure. At this point, alternative characterization (such as DMA, with the NETZSCH DMA 242 *Artemis*) would be more suitable; the disposable plate system is designed to monitor the progression of the cure.



5 Example curing data for two different samples

<sup>1</sup> There are multiple definitions of the gel point. A simple convenient one (as detailed in ASTM 4473) is where G' = G'' (i.e., phase angle = 45°) but this would change depending on the applied frequency of the oscillation. A more detailed definition of the sol-gel transition is detailed by Winter & Chambon (Chambon, 1987) as the point where the phase angle is independent of frequency (rather than 45°).



Use of the Kinexus Disposable Plate Systems

During curing processes, some materials can change in volume, either expansion of contraction, which can also be monitored on Kinexus during the same curing test. Often, these changes are unwanted. However, sometimes (depending on end use) this can provide helpful information on the process, such as sample expanding in a cavity that the adhesive is required to fill, or avoid contraction which can put immediate and unwanted stresses onto a bond.

For any of these scenarios, Kinexus is available to quantity this change. Either by fixing the gap and monitor the changes in "Normal Force" (axial force; positive value being up, negative down) to indicate the "tendency" of a material to want to change volume, or maintaining.

# Advanced Test Control: rSpace Triggers

The Kinexus rSpace software has advanced control over normal force measurements by use of the "Trigger" system. This means, even during a measurement, the control of the sequence can change with the trigger condition met, hence "triggering" a different action, for example.

Figure 6 shows a typical example of how the trigger system works for a curing material which is initially "viscously dominated" (i.e., a viscoelastic liquid, so a tendency to flow). As such, an initial applied normal force would push (flow) the sample out of the measurement gap. In this viscous state, the sample will flow as it cannot apply an outwardly force as the sample will just otherwise deform (to relax) instead (hence no net change in normal force). As such, the sequence starts in a fixed gap measurement, "triggering" to normal force control (i.e., variable gap) once the sample is too high in modulus so no longer likely to be viscously dominated enabling the sample to likely to expand or contract.

Further bespoke and flexible control, is available with this unique and flexible trigger system supplied with rSpace; most instrument variables are available to be trigger on.

The Kinexus disposable plate system retains compatibility for use with the Kinexus solvent trap systems:

- Kinexus passive solvent trap cover for standard plate and cylinder cartridges only, made of stainles steel
- Kinexus active solvent trap cover for active hood cartridge only, made of stainless steel

### References

Chambon, F. A. (1987). Linear viscoelasticity at the gel point of a crosslinking PDMS with imbalanced stoichiometry. Journal of Rheology 31.8, 683-697.



6 Trigger control to change measurement conditions during a sequence

