


APPLICATION NOTE

Pharmaceuticals – TG 209 F1 Libra-FT-IR

Improving Process Conditions – Determining the Thermal Stability of Eudragit® by Means of TGA-FT-IR



1 Structure of the monomer of Eudragit® L100-55 [1]

Introduction

Eudragit® is the trade name of polymethacrylate-based co-polymers used to target the release of a drug in the desired parts of the gastro-intestinal tract. Eudragit® exists in several compositions that differ from each other in terms of the functional groups located on the side chains.

This results in differing dissolution behavior depending on the pH value of the environment. For example, Eudragit® L100-55 (figure 1) is soluble in intestinal fluids from pH 5.5 upward, but is not soluble in gastric fluids with lower pH values. Thus, it is used in combination with drugs to be released in the duodenum after passing the stomach [1, 2, 3].

Thermal analysis of Eudragit® products is crucial for different reasons: Firstly, they differ from each other in their glass transition temperature (T_g). Even Eudragit® polymers with a similar chemical composition show differences in T_g , depending on the monomer ratios [1]. Thus, the determination of T_g allows for identification

of the different Eudragit® polymers. Secondly, optimum process conditions, e.g., for hot melt extrusion, require knowledge of the glass transition temperature and thermal stability of the polymer [3].

For this reason, the decomposition process of Eudragit® L100-55 (Evonik Industries) is investigated by means of a thermobalance (TGA) coupled to an FT-IR spectrometer.

Measurement Conditions

The TGA-FT-IR measurement was performed using a NETZSCH TG 209 **F1** Libra thermobalance. To investigate and identify the gases released during thermogravimetric analysis, they were transferred directly into the gas cell of an FT-IR system by Bruker Optics.

The measurement was carried out on 7.33 mg of Eudragit® L100-55, using an open aluminum oxide crucible. The sample was heated between room temperature and 600°C at 10 K/min in a nitrogen atmosphere (40 ml/min).

Measurement Results


Figure 2 depicts the mass changes of Eudragit® L100-55 between 40°C and 600°C. The first mass-loss step of 0.8% indicates the release of surface water up to 100°C. The second mass loss of 5.9% at 200°C (DTG peak) is also associated with the release of water confirmed by the FT-IR spectrum (figure 3). The temperature of the process indicates the release of crystal water. In addition, bands occur in the wavelength range 3000 - 2800 cm^{-1} and above 1000 cm^{-1} . These bands represent CH_2 and CH_3 molecules that indicate the start of decomposition of the Eudragit® sample.

APPLICATIONNOTE Improving Process Conditions – Determining the Thermal Stability of Eudragit® by Means of TGA-FT-IR


The peak at 294°C in the DTG curve is associated with yet another step in the decomposition process: the release of carbon dioxide and probably ethanol (figures 4 and 5). This can be explained by the splitting of an ester group off the Eudragit® molecule (figure 6).

The last and main decomposition step, with a mass loss of 88.5%, occurs at 393°C (DTG peak temperature). The characteristic bands of ethanol and of carbon dioxide can still be detected in the FT-IR-spectrum of the gases

released at 393°C (figure 7). In addition, carbon monoxide (2300 cm^{-1} to 2100 cm^{-1}) and an ester substance are present, which can be seen in the carbonyl band at 1749 cm^{-1} . It suggests that the ester part $\text{C}_2\text{H}_5\text{-O-CO-C}_x\text{H}_y$ breaks off from the molecule (see red indication in figure 10). The two vibration bands at 1460 cm^{-1} and 1380 cm^{-1} are probably due to parts of the carbon backbone. As an example, a comparison spectrum of ethyl acetate and 3-methyloctane is illustrated in figures 8 and 9.




2 Mass changes of Eudragit® L100-55 during heating to 600°C




3 FT-IR spectrum of the gases released at 206°C during heating of Eudragit® L100-55

APPLICATIONNOTE Improving Process Conditions – Determining the Thermal Stability of Eudragit® by Means of TGA-FT-IR




4 Comparison of the FT-IR spectrum of the gases released at 295°C (top) with the EPA-NIST-FT-IR spectrum of carbon dioxide (bottom)




5 Comparison of the FT-IR spectrum of the gases released at 295°C (top) with the FT-IR spectrum of ethanol (bottom [4])

APPLICATIONNOTE Improving Process Conditions – Determining the Thermal Stability of Eudragit® by Means of TGA-FT-IR




- 6 The split-off of an ester could explain the detection of ethanol in the gases released at 295°C




- 7 FT-IR spectrum of the gases released at 393°C

APPLICATIONNOTE Improving of Process Conditions –Determining the Thermal Stability of Eudragit® by Means of TGA-FT-IR




8 Comparison of the FT-IR spectrum of the gases released at 393°C (top) with the EPA-NIST-FT-IR spectrum of ethyl acetate (bottom)



9 Comparison of the FT-IR spectrum of the gases released at 393°C (top) with the FT-IR spectrum of 3-methyloctane (bottom) [4]

APPLICATIONNOTE Improving of Process Conditions – Determining the Thermal Stability of Eudragit® by Means of TGA-FT-IR



10 The Break-up of the ester group

Conclusion

The start of decomposition of Eudragit® is closely related to the thermal stability. It results in changes in sample mass during storage or thermal treatment. The mass changes can be identified by means of thermogravimetry. However, clear identification of the released gases – and thus a reliable interpretation of the mass losses – is only possible when the thermobalance is coupled to an FT-IR device. This allows for reliable conclusions to be drawn as to whether a given mass loss can be attributed to decomposition or merely to the release of water.

Under the selected conditions (inert atmosphere, heating rate of 10 K/min), the investigated Eudragit® sample starts to decompose at 185°C (onset temperature of the TGA curve). The fact that this is really the start of decomposition is revealed by the occurrence of C-H bonds in addition to crystal water.

Literature

- [1] Investigation of thermal and viscoelastic properties of polymers relevant to hot melt extrusion – III: Polymethacrylates and polymethacrylic acid based polymers.
Tapan Parikh, Simerdeep Singh Gupta, Anuprabha Meena, Abu T.M. Serajuddin
- [2] Eudragit and its Pharmaceutical Significance, Satish Singh Kadian, S.L. Harikumar
- [3] http://healthcare.evonik.com/sites/lists/nc/documentshc/evonik-eudragit_brochure.pdf → [LINK](#)
- [4] <https://webbook.nist.gov/>