

## Smart Thermal Analysis: Measurements Wanted?

Dr. Alexander Schindler



### Introduction

In day-to-day laboratory work, there are typical questions that often arise. For a new sample regarding analysis, what are the suitable measurement conditions such as temperature program, sample mass or the right crucible? And what measurement results can be expected? Perhaps such a sample was already measured by you in the past – or maybe by NETZSCH.


Wouldn't it help a lot to simply search in a database for thermal analysis? *Identify*, which is a part of the *Proteus*® analysis software, is the solution!

The primary function of the *Identify* database system is the automatic recognition and comparison of measurement curves. This serves, for example, for quality control and failure analysis. Pure data mining (storing, searching and finding of data) is, of course, the second main application.

**How to Carry Out Data Mining with *Identify*?**

One approach in using *Identify* is to search for similar database measurements or literature data on the basis of an input measurement curve. This is illustrated in figure 1, where a Thermomechanical Analysis (TMA) measurement is compared with various database curves for the same material. The similar TMA curve was found automatically by *Identify*, while the DSC and TGA curves in

this case were found simply by alphabetical sorting of all database measurements shown. It is interesting to see that the glass transition occurs between -70° and -60°C (see also the Differential Scanning Calorimetry curve) and that the decomposition – which should normally be avoided in DSC or TMA measurements – starts slowly at above 150°C; this can be seen in the Thermogravimetry (TGA) curve.



1 *Identify* database search based on a TMA curve of an NR (natural rubber) sample. Overlaid are a similar TMA curve, a DSC and a TGA curve of NR found in the database.

## APPLICATIONNOTE Smart Thermal Analysis: Measurements Wanted?

The second approach in data mining – in which no similar input measurement is required – is related to the “Manage Libraries/Classes” function (see figure 2). In this case, “NR” was typed manually into the search field for measurements and literature data, and consequently,

the three different NR measurements were found. Right-clicking restores the database measurement to also reveal details of the measurement conditions such as temperature program, sample mass, purge gases and sample crucible.

**Libraries**   Folders   Classes

Libraries

Name	Entries
Pharma_Food_Cosmetics NETZSCH	239
Polymers DSC KIMW	1100
Polymers NETZSCH	190
Polymers Poster NETZSCH	70

Measurements/Literature Data

NR

Name	Effects	DIL	Cp
NR_DIL	0	x	
NR_DSC	?		
NR_TGA			

Effects

Name
Glass Transition at -61,71°C (Mid)
Endothermic Effect at 24,17°C (Onset)

DIL/TMA and Cp Data

Name
------

Close

2 Searching within the *Identify* database using the name filter (indicated by the red circle)


## APPLICATIONNOTE Smart Thermal Analysis: Measurements Wanted?

Last but not least, *Identify* also offers the ability to filter according to the heating rate and the sample mass, as well as according to whether glass transitions or endothermic and exothermic effects were evaluated. It is also possible to filter according to the occurrence, or lack of occurrence, of a string of letters in a measurement name (see figure 3).

fields (ceramics, inorganics, metals, alloys, organics, pharma, food, cosmetics and polymers). Available as an option is the [KIMW database](#) developed by the Kunststoffinstitut Lüdenscheid, Germany, with DSC curves for 1000 different commercially available polymer grades; here information about the polymer supplier, color and filler material/content is also available.

### The Database Content of *Identify*

Shown in figure 4 are the NETZSCH libraries, currently containing 1294 entries, which cover various application



3 Filter settings within *Identify*

Library	Entries
+ Literature Data Poster NETZSCH	248
Ceramics_Inorganics NETZSCH	302
Metals_Alloys NETZSCH	143
Organics NETZSCH	172
Pharma_Food_Cosmetics NETZSCH	239
Phase Change Materials	14
Polymers DSC KIMW	1000
Polymers NETZSCH	190

NETZSCH

User

KIMW

4 *Identify* database contents: NETZSCH libraries, which are always included (1294 entries), the optional KIMW library (1000 entries), and in an example of a library created by the user.

## APPLICATIONNOTE Smart Thermal Analysis: Measurements Wanted?

### Summary

Any measurement stored in *Identify* can be easily found and its associated measurement conditions and evaluations retrieved. This treasure chest of information can be helpful in advance of, or following, any measurement.

### References

This article is also available as a blog:

<https://ta-netzsch.com/smart-thermal-analysis-measurements-wanted>

These previous released blog articles might also be of interest:

[Smart Thermal Analysis \(Part I\): AutoEvaluation](#) of DSC, TGA and STA curves

[Smart Thermal Analysis \(Part II\)](#): Identification of Measurements via Database Search

[Smart Thermal Analysis \(Part IIb\)](#): *Identify* ... the Most Comprehensive Database in Thermal Analysis

[Smart Thermal Analysis \(Part III\)](#): *AutoEvaluation* of DIL and TMA curves