

Technical Specifications

TDW 4040

Measuring range	R: 0.10 to 8.00 m ² ·K/W
Specimen size (L x W)	1500 mm x 1500 mm (brick wall)
Specimen thickness (H)	Up to 360 mm
Temperature range	<ul style="list-style-type: none">■ Cold section: -10°C to 40°C■ Warm section: +10°C to 40°C
Interface	1x RS 232, 1x Gigabit Ethernet
Dimensions (H x W x D)	300 cm x 500 cm x 360 cm
Power supply	230 V / 400 V / 50 Hz
Weight	3890 kg